


# MicroPak® Classic

## Protecting Your Investment

The MicroPak® Classic Series provides a simple and convenient solution to the organization, sterilization, and transportation of microsurgical instrument sets while providing a high level of protection.

### Benefits


- Provides maximum protection as well as ease of fixation for a variety of instruments.
- Withstands repeated autoclave and dry heat to 320 F (160 C) as well as EO.



### MicroPak Classic Series I, Base, Lid and Mat

Description	Part #
Base, Lid and Mat, 9.6 x 5.5 x .7 inches (24.4 x 14 x 1.8 cm)	0105
Base, Lid and Mat, 12.8 x 7.3 x .7 inches (32.5 x 18.5 x 1.8 cm)	0116
Base, Lid and Mat, 18 x 10.3 x .7 inches (45.7 x 26.2 x 1.8 cm)	0119


*\*All dimensions are internal*



### MicroPak Base, Lid and Mat

Description	Part #
Base, Lid and Mat, 5.5 x 9.5 x 2.5 inches (14 x 24.1 x 6.4 cm)	1015


*\*All dimensions are internal*



## MicroPak Classic Series II, Base, Lid, Insert Tray and Mat

Description	Part #
<b>Base, Lid, Insert Tray and Insert Mat, 8.5 x 4.5 x 1.5 inches (21.6 x 11.4 x 3.8 cm)</b>	<b>0207</b>
<b>Base, Lid, Insert Tray and Insert Mat, 9.5 x 8 x 1.5 inches (24.1 x 20.3 x 3.8 cm)</b>	<b>0212</b>
<b>Base, Lid, Insert Tray and Insert Mat, 11.8 x 6.3 x 1.5 inches (30 x 16 x 3.8 cm)</b>	<b>0217</b>
<b>Base, Lid, Insert Tray and Insert Mat, 17.8 x 9.7 x 1.5 inches (45.2 x 24.6 x 3.8 cm)</b>	<b>0220</b>

*\*All dimensions are internal*



**SYMMETRY Medical<sup>®</sup>** Inc.

27 Wrights Landing, Auburn, Maine 04210-9915, USA  
T: +1.800.245.3300 +1.207.786.2775 Fax: +1.207.786.5553

Grand-Rue 53, Coregemont, Switzerland, CH 2606  
T: +011.41.32.488.1010 Fax: +011.41.32.488.1000  
[www.symmetrymedical.com](http://www.symmetrymedical.com)