

Acromioplasty Retractor

The Acromioplasty Retractor is used to retract and protect the humeral head during resection of the inferior acromial surface.

The two prongs hook the posterior aspect of the acromion for retraction while the file is used to smooth rough edges on the acromion post-resection.

BENEFITS

- Retractor protects the humeral head during acromioplasty.
- Dual pronged working end improves exposure of acromion via secure retraction of the humeral head.
- File end allows smoothing of the acromion after resection.

Contact your local Sales Manager for more information.



DESCRIPTION	PART #
Acromioplasty Retractor	13-008-000



SYMMETRY Medical[®] Inc.

486 West 350 North, Warsaw, IN 46582, USA
T: +1.574.267.8700 Fax: +1.574.267.7306
www.symmetrymedical.com